

■ Light 2 Unisepta®



Ready to use, alcohol based solution for the rapid disinfection of surfaces of medical devices.
With an action time from just 1 minute, Unisepta® Light 2 is ideal for use between patients.

Key product features

- Ethanol based
- Effective from 1 minute
- Leaves no residue
- Fresh and light perfume
- CE marked – Complies with Medical Devices Directive 93/42/EEC
- Only 40% alcohol



USF Healthcare S.A. | Rue François Perréard 4 | CH-1225 Chêne Bourg | Geneva | Switzerland
Tél + 41 22 839 79 00 | Fax + 41 22 839 79 10 | info@usfhc.com
www.usfhealthcare.com



Light 2 Unisepta®


Protocol for use

Before use


Ensure appropriate gloves are worn.



Ensure that the area to be treated has been previously cleaned and is resistant to alcohol.



Spray the area directly or in to a non-woven cloth, then wipe. Ensure that the entire area is treated and comes in to contact with the solution.




Respect the contact time indicated for the desired antimicrobial activity.

Do not rinse!

Precautions

Dangerous: respect the precautions for use (established according to the European regulations in force for classification and labelling of chemicals). Before each use, read the label and the information concerning the product. Storage between +5°C and +25°C. Reserved for professional use. Medical device Class IIa (complies with medical devices directive 93/42/EEC modified).

Efficacy



* PRV: Surrogate of Hepatitis B

* BVDV: Surrogate of Hepatitis C

The Robert Koch Institute recommends the use of disinfectants claiming limited virucidal activity to interrupt the infection-chain. According to the RKI, limited virucidal activity (effective against enveloped viruses) corresponds to an activity against Vaccinia virus and BVDV.

Conforms to European standards

Bactericide: EN 1040, EN 13727, EN 13697

MRSA: EN 13727

Yeasticide: EN 1275, EN 13624, EN 13697

Active on viruses according to EN 14476

Product composition

100g of solution contains: Ethanol (327 mg/g), N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine (1.4 mg/g), didecyldimethylammonium chloride (1.1 mg/g), perfume and excipients.

Physical data

Colour: Colourless

Perfume: Light perfume

Shelf Life: 2 years

Container: Made of recyclable HDPE material

Packaging

750ml Spray	1L	2,5L	5L
✓			✓



ECOEngage
Our Ecological Engagement



USF Healthcare
Unident • Steridine • Farmec