

# Steridine® Ultra Clean

Enzymatic detergent for cleaning of surgical and medical instruments, thermosensitive instruments and endoscopic equipment.



- Aldehyde free
- Economical dilution of 0.5%
- Protease enzymatic formulation for optimum cleaning
- Cares for delicate and expensive instrumentation
- Rapid cleaning action - effective in just 1 minute
- Fresh citrus fragrance
- Not classified dangerous for the environment when diluted

## Product composition

Surfactants, sequestering agent, protease, stabilising agent, perfume, excipients.

## Physical data

Colour: Green  
Perfume: Citrus  
pH: 7.6 when diluted at 0.5% / 8.5 undiluted  
Container: Made of recyclable HDPE material

## Packaging

12 x 1L ..... ref. TE00110  
2 x 5L with 1 x 25ml pump ..... ref. TE00114

## Precautions

Dangerous – respect the precautions for use (established according to the European regulations in force for classification and labelling of chemicals). Before each use, read the label and the information concerning the product. Storage between +5°C and +35°C. Reserved for professional use. Medical device Class I (complies with medical devices directive 93/42/EEC modified).