

# Neutraliser Micro 10<sup>®</sup> WD



Rinsing and neutralising liquid, with lubricating effect, for dental and medical instrumentation processed in automatic washer disinfectors. Double action to ensure quick drying without trace.

## Key product features

- Excellent material compatibility
- Lubricating effect
- Can be used with all grades of rinsing water, treated or untreated
- Compatible with all types of washer disinfectors
- Injection diluted between 0.03% to 0.3% depending on type of machine
- CE marked – Complies with Medical Devices Directive 93/42/EEC
- Use in combination with Micro 10<sup>®</sup> WD Detergent



USF Healthcare S.A. | Rue François Perréard 18 | CH-1225 Chêne Bourg | Geneva | Switzerland  
Tél + 41 22 839 79 00 | Fax + 41 22 839 79 10 | [info@usfhc.com](mailto:info@usfhc.com)  
[www.usfhealthcare.com](http://www.usfhealthcare.com)



# Micro Neutraliser 10<sup>®</sup> WD

## Product composition

Non ionic and anionic surfactants, isopropanol, excipients.

## Physical data

Colour: Clear

pH: 2 to 3 (pure product at 20°C)  
Approximately 7 (diluted at 0.2%)

Shelf life: 3 Years

Container: Made of recyclable HDPE material

## Packaging

| 1L | 2,5L | 5L |
|----|------|----|
|    |      | ✓  |

## Precautions

Dangerous: respect the precautions for use (established according to the European regulations in force for classification and labelling of chemicals). Before each use, read the label and the information concerning the product. Storage between +5°C and +25°C. Reserved for professional use. Medical device Class I (complies with medical devices directive 93/42/EEC modified).

