



Septoscrub

Antiseptic liquid soap for **EXTERNAL USE**

Technical Information

1. Composition

100g of solution contains:

Active ingredient:

IRGASAN DP 300 (2,4,4'-trichloro-2'-hydroxydiphenyl ether – CAS No. 3380-34-5) 0.5g

Excipients:

Surfactants, essence, preservative, coformulants and purified water up to 100g

2. Product characteristics

SEPTOSCRUB is a ready-to-use aqueous solution with antiseptic and detergent properties. It has a pH of 5.5, that is perfectly compatible with the skin.

The product has a broad spectrum of biocidal action, which also includes viruses (EN 14476).

The efficacy of SEPTOSCRUB has been tested for hygienic handwash, as stated by EN 1499 European standard.

Mechanism of action:

The active ingredient of SEPTOSCRUB reduces surface tension at the point of absorption and is rapidly absorbed on surfaces with a negative charge (proteins, bacteria etc.). It has precipitating, complexing and denaturing effects on proteins and causes enzymatic modifications to the cytoplasmic membrane of micro-organisms.

3. Range of application and instructions for use

3.1 Scope of use

- Hygienic handwash.
- Complete washing of patient's body before surgery.

3.2 Instructions for use

SEPTOSCRUB is a ready-to-use antiseptic liquid soap:

For hand antiseptics:

- Wash the hands under running water;
- Apply approximately 7 ml (2 pumps of the dispenser) of the solution onto the palm of your hand;
- Rub the solution systematically over the hands and wrists for around 1 minute;
- Lastly, rinse under running water.

Since the active ingredient IRGASAN DP 300 is incompatible with chlorinating substances (see paragraph 3.6), it is recommended that you use tap water with low chlorine content for wetting your hands.

3.3 Special precautions for use

THE PRODUCT IS FOR EXCLUSIVE EXTERNAL USE.

The use of the product, even if prolonged, has not lead to hypersensitivity.

3.4 Use during pregnancy and breastfeeding

There are no limits on use during pregnancy or breastfeeding.

3.5 Interactions with other products

Avoid using at the same time other antiseptics and/or detergents.

3.6 Incompatibility

The substances in the formula are incompatible with cationic detergents and chlorinating compounds.

4 Packaging

Bottle

20 x 500ml ref. PF04987

5 Quality Control

The components (raw materials, containers, labels etc.) and production steps of each individual batch are carefully controlled in compliance with the SOP (Standard Operative Procedures) required by standard UNI EN ISO 9001 and Good Manufacturing Practices.

6 Storage and stability

Keep the container properly closed, away from light and heat sources.

SEPTOSCRUB in the original sealed packaging has a validity of 36 months.

If the package is correctly opened and closed at the end of each disinfection operation, without contamination of the residual part, the product maintains its chemical-physical and microbiological characteristics for 90 days.

7 Regulatory references

Holder of the registration:

Nuova Farmec S.r.l. – Via Flemming, 7
37026 Settimo di Pescantina (VR) – Italy.

Presidio Medico Chirurgico authorized by Italian Ministry of Health registration N. 16669.



INFORMATION RESERVED FOR HEALTH CARE PROFESSIONALS AND PROFESSIONAL USERS
